

FOR TRUSTWORTHY BP MONITORING AT HOME

TRUSTCHECKTM BPM 3.0

AUTOMATIC UPPER ARM BLOOD PRESSURE METER



- WHO INDICATOR
- INTELLIGENT PRESSURIZATION
- BIG LCD SCREEN
- 99 TESTS MEMORY FOR 2 USERS EACH
- IRREGULAR HEARTBEAT DETECTION
- FACILITY FOR MICRO USB POWER ADAPTER



TRUSTCHECK™ BPM 3.0

AUTOMATIC UPPER ARM BLOOD PRESSURE METER



WHO INDICATOR

Makes it easy to understand Results



INTELLIGENT Pressurization

Oscillometric Method offers pressurization as per intensity of Blood Pressure



BIG LCD Screen 7.6 x 5.2 cm

For ease of reading the results



99 MEMORIES for 2 Users

Store & refer to userwise results easily



IRREGULAR HEARTBEAT DETECTION

Can help in detection of pulse irregularities at an early stage



Facility for MICRO USB POWER Adapter

Safe & Economical




Stages of Blood Pressure levels		Systolic (in mmHg)	Diastolic (in mmHg)
Grade 3	Severe hypertension	≥180	≥110
Grade 2	Moderate hypertension	160 ~ 179	100 ~ 109
Grade 1	Mild hypertension	140 ~ 159	90 ~ 99
High-Normal		130 ~ 139	85 ~ 89
Normal		120 ~ 129	80 ~ 84
Optimal		<120	< 80

*Source: WHO,2003

BOX CONTENTS:

- 1 Blood Pressure Meter
- 1 Cuff
- 4 x AAA Batteries
- 1 Carry Pouch
- 1 User Manual

For more information contact:



TECHNICAL SPECIFICATION

Product Type	Automatic Upper arm type: ARK-BP-04
Measurement Range	DIA: 40-199mmHg; SYS: 60-260mmHg Pulse 40-180beats/min
Accuracy	Blood pressure SYS/DIA : ±3mmHg Pulse: < 60: ± 3%; >60(incl.): ±5% of the reading
Memory	2 x 99 groups memory (2 users)
Measuring Method	Oscillometric method
Device Weight	270 grams
Cuff Circumference	22 - 32 cm
Power Source	Internal-DC 6V (4*AAA) / External-DC 5V 1A
Operating Conditions	Temperature: 10°C to 40°C Humidity: < 85%RH
Storage And Shipping Conditions	Temperature: -20°C to 50°C Humidity: < 85% RH

Follow us on Social media for interesting info & updates



Diabetes Management
Mobile App Scan QR Code
to Download (IOS & Android)



ARKRAY



arkray_official



Customer Support (Mon-Fri:9:30 am to 6:30 pm)

WhatsApp: +91 77770 53311 | Email: customersupport@arkray.co.in

ARKRAY Healthcare Pvt. Ltd | 701-702, Opulence, 6th Road, TPS III, Santacruz East, Mumbai 400055, Maharashtra, India